



JENWAY

Spectrophotometer
Accessories



Jenway® Spectrophotometer Accessories

Jenway® spectrophotometers are very versatile and flexible instruments with an extensive range of accessories available which have been specifically designed to be easily interchangeable. The active accessories include the 8-cell automatic turret, 6 cell water heated automatic turret, sipper pump, Peltier and combined sipper/Peltier. Jenway also offer a range of passive accessories that can suit any application or budget.



660201

660501

660401

660301

660701



685131

685304

685204

685005

685401

Ordering information

Model number	Description	For use with
660201	Sipper pump	67 series
660301	Electronic Peltier for temperature controlled applications	67 series
660401	Automatic 8 position cuvette holder	67 series
660501	Water heated automatic six-cell changer	67 series
660701	Combined Sipper/Peltier	67 Series
661201	Single water heated 10x10mm cell holder	67 Series
685131	Water heated 10x10mm single cell holder	6850
685304	Micro-cuvette holder	6850
685204	10x10mm path length cuvette holder	6850
685005	10 to 100mm cuvette holder	6850
685401	8 position automatic cell changer	6850

630204



637071



630005



630304




736201



735401




735201



735301



735701



660101



543001



Ordering information

Model number	Description	For use with
630204	10x10mm cuvette holder	63, 73 and 67 series
637071	16/24mm test tube holder	63, 73 and 67 series
630005	Variable path length cuvette holder	63, 73 and 67 series
630304	Micro-cuvette holder with reduced aperture	63, 73 and 67 series
634001	Manual 4 position cuvette holder	63 Series
648001	Water heated single cuvette holder	63 Series
543001	40 column serial printer	63 Series
736201	Water heated 10 x 10mm path length cuvette holder	73 Series
735401	Automatic 8 position cuvette holder	73 series
735201	Sipper pump	73 series
735301	Electronic Peltier for temperature controlled applications	73 Series
735701	Combined Sipper/Peltier	73 series
660101	Internal Printer	67 and 73 Series

Cuvettes



Jenway offer a range of glass and quartz cuvettes as well as disposable plastic cuvettes. The glass and quartz cuvettes are available in a range of path lengths from 10 to 100mm. Jenway also offer a range of glass and quartz flow-through and micro-volume cuvettes. Glass cuvettes are suitable for measurements in the visible region whereas quartz cuvettes may be used in both the UV region and visible region.

Part Code	Material	Description	Visible wavelength	UV/Visible wavelengths	Path length	Fill volume
060 084	Plastic	Pack of 100 plastic cuvettes	✓	-	10mm	2.4 to 4.5ml
060 229	Plastic	Pack of 500 plastic cuvettes	✓	-	10mm	2.4 to 4.5ml
060 087	Plastic	Pack of 100 plastic cuvettes	✓	-	10mm	1.5 to 3.0ml
060 230	Plastic	Pack of 100 plastic cuvettes	✓	✓	10mm	2.5 to 4.5ml
060 232	Plastic	Pack of 500 plastic cuvettes	✓	✓	10mm	2.5 to 4.5ml
060 231	Plastic	Pack of 100 UV plastic cuvettes	✓	✓	10mm	1.5 to 3.0ml
035 143	Plastic	Pack of 100 plastic cuvettes, UV	✓	✓	10mm	70µl to 1.5ml
035 027	Glass	10mm path length cuvette	✓	-	10mm	2.0 to 3.5ml
035 086	Glass	20mm path length cuvette	✓	-	20mm	Up to 7.0ml
035 029	Glass	40mm path length cuvette	✓	-	40mm	Up to 14.0ml
035 087	Glass	50mm path length cuvette	✓	-	50mm	Up to 17.5ml
035 079	Glass	100mm path length cuvette	✓	-	100mm	Up to 35.0ml
035 123	Glass	Micro (500µl) cuvette	✓	-	10mm	500µl
035 126	Glass	Semi-micro (1ml) cuvette	✓	-	10mm	1ml
035 045	Glass	Flow-through cell	✓	-	10mm	80µl
035 044	Quartz	Flow-through cell	✓	✓	10mm	1.8ml
035 047	Quartz	Flow-through cell	✓	✓	10mm	80µl
035 028	Quartz	10mm path length cuvette	✓	✓	10mm	2.0 to 3.5ml
035 056	Quartz	20mm path length cuvette	✓	✓	20mm	Up to 7.0ml
035 030	Quartz	40mm path length cuvette	✓	✓	40mm	Up to 14.0ml
035 281	Quartz	50mm path length cuvette	✓	✓	50mm	Up to 17.5ml
035 282	Quartz	100mm path length cuvette	✓	✓	100mm	Up to 35.0ml
035 124	Quartz	UV ultra micro (50µl) cuvette	✓	✓	10mm	50µl
035 125	Quartz	UV micro (500µl) cuvette	✓	✓	10mm	500µl
035 127	Quartz	UV semi-micro (1ml) cuvette	✓	✓	10mm	1ml
035 138	Quartz	UV ultra-micro (100µl) cuvette	✓	✓	10mm	100µl

Data management software

Bibby Scientific's new data management software offers a powerful combination of functional but simple applications to store and retrieve laboratory test results and information. Qualoupe Lite is an ideal data storage solution for research groups, commercial laboratories and small to medium enterprises. Qualoupe Lite has an in-built option to change the displayed text to English, French or German language.



Product Code	Description
QUALITE	Qualoupe Lite Data Management Software

JENWAY
A Bibby Scientific Brand

Bibby Scientific - UK (Group HQ)
Beacon Road,
Stone, Staffordshire,
ST15 0SA,
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
e-mail: sales@bibby-scientific.com

sales@bibby-scientific.com

www.jenway.com



Find out more!!!

Please scan the QR/Mobile Tag with your smartphone for more information